

---

---

## **Cable-Ready™ Cable Grip System Bone Plate**

JACK JENNINGS, BSME, LYNN KIRKPATRICK, BSME, PE, AVIJIT MUKHERJEE, MSME, MBA\*

*Investigation performed at the Research Laboratories of Zimmer, Inc., Warsaw, Indiana,  
and at Pioneer Laboratories, Marquette, Michigan*

*The Cable-Ready Cable Grip System was developed by Pioneer Laboratories, Marquette, Michigan*

---

---

**Summary:** Shaft fractures near the prosthesis in patients who have previously received a total hip, knee, or shoulder prosthesis are serious complications. Management of these fractures can present a unique challenge for the orthopaedic surgeon. Open reduction and internal fixation with plates and screws and/or cerclage wires or cables is a well established treatment modality. The use of Ogden plates and Parham bands validated these techniques many years ago. The *Cable-Ready* Cable Grip System bone plate is an evolution of previous fixation concepts providing immediate rigid fixation of periprosthetic long bone fractures and allowing early mobilization of the affected limb. This paper provides information and comparative test data for the *Cable-Ready* bone plate, Zimmer *ECT*® Internal Fracture Fixation Broad Compression Plate, and other competitive devices (such as Howmedica, Dall-Miles, and Biomet).

### **Uses and Indications**

The *Cable-Ready* bone plate is indicated for use where wire, cable, or band cerclage is used in combination with bone plates to provide fixation and/or stabilization of long bones. Other uses include fracture below or above a prosthesis, comminuted femoral shaft fractures, nonunion of fractures with failed hardware, reconstruction of bone defects, and total hip or knee revisions with bone loss or fracture.


### **Design**

The *Cable-Ready* bone plate was designed and developed to address complicated fractures of the femur, tibia, and humerus. The plate design makes the cables an integral part of the plate. No additional external crimping devices are required. Each cable is passed through the bone plate and the cable is locked by turning a set screw which crimps

\*Jack Jennings is a Research Engineer working in the Fatigue and Fracture Mechanics Laboratory at Zimmer. He holds a B.S. Degree in Materials Engineering from Purdue University. He is a member of the American Society of Metals, Society of Experimental Mechanics, and Society of Advanced Materials and process Engineering.

Lynn Kirkpatrick is a senior Research Engineer, Zimmer, where he is responsible for mechanical testing in the support of the research and development and production functions. He earned his B.S. Degree in Mechanical Engineering from Tri-State University in 1974 and is a licensed Professional Engineer (PE) in Indiana. Mr. Kirkpatrick has 15 years of experience in the orthopaedic industry.

Avijit Mukherjee is a group manager in Fracture Development area of the Orthopaedic Implant Division of Zimmer. He holds a M.S. degree in Mechanical Engineering from Tennessee Technological University and an M.B.A. degree from the University of Missouri, St. Louis. Mr. Mukherjee has over 20 years of experience in materials research, product development, and program management in defense, aerospace, and orthopaedic industries.



the cable securely in place. This provides a very low profile device. The fixation of the bone plate is very effective at resisting biomechanical loads in flexion and extension. Bone screws can be inserted through the plate to resist torsional loads. Additional bone plate information including materials and sizes is shown in Table 1.

Product testing was conducted on the cable plate and other similar devices including competitive plates. The unique design features of the *Cable-Ready* bone plate have contributed to significantly higher static bending strength and stiffness compared to the Howmedica, Dall-Miles, and *ECT* Internal Fracture Fixation bone plates. The four-point bend test results were performed according to ASTM guidelines and are summarized in Table 2. The Zimmer bone plate and *ECT* plate were also tested in fatigue (repeated cyclic loading). The results of the fatigue test are summarized in

TABLE 1  
CABLE-READY CABLE GRIP BONE PLATE  
MATERIALS AND SIZES

Item	Material	Size
Bone Plate	316L stainless steel	6-hole: 7.36" (187mm) 8-hole: 9.69" (246mm) 10-hole: 12.01" (305mm)
Cable	316L stainless steel	1.8mm diameters X 24" (610mm) long
Set Screw	22-13-5 stainless steel	3.8mm long X 3.0mm hex
Tensioning Bit	17-4 stainless steel	

TABLE 2  
4-POINT STATIC BEND TEST RESULTS  
(ASTM F-382)

Bone Plate Type	Bone Plate Material	Bending Stiffness	Bending Strength
Zimmer <i>Cable Ready</i>	316L SST	10,432 lb-in/in	834 lb-in
Howmedica Dall-Miles Cable Plate	316L SST	8410 lb-in/in	483 lb-in
Zimmer 9-hole <i>ECT</i> Broad Compression Plate	316L SST	6063 lb-in/in	561 lb-in

Table 3. In addition, plate profiles (thicknesses) were measured for the Zimmer bone plate, Dall-Miles, and Biomet bone plates. These results are shown in Table 4.

### Conclusion

- The average bending stiffness of *Cable-Ready* bone plates is 72 percent greater than 9-hole *ECT* broad compression plates and 24 percent greater than Dall-Miles plates.
- The average bending strength of *Cable-Ready* bone plates is 49 percent greater than *ECT* broad compression plates and 73 percent greater than Dall-Miles plates.

TABLE 3  
4-POINT FATIGUE BEND TEST  
(@ 30HZ, 2 MILLION CYCLES)

Bone Plate Type	Bone Plate Material	Fatigue Endurance Stiffness
Zimmer <i>Cable Ready</i>	316L SST	
Zimmer 9-hole <i>ECT</i> Broad Compression Plate	316L SST	140 lb

### Conclusion

- Fatigue endurance strength limit of *Cable-Ready* bone plates is 43 percent higher than 9-hole *ECT* Broad Compression plates.

TABLE 4  
4-POINT FATIGUE BEND TEST  
(@ 30HZ, 2 MILLION CYCLES)

Bone Plate Type	Bone Plate Material	Cable/Plate Profile
Zimmer <i>Cable Ready</i>	316L SST	7.5mm
Howmedica Dall-Miles Cable Plate	316L SST	9.1mm
Biomet Cable Plate	316L SST	

### Conclusion

- Profile of *Cable-Ready* bone plates is lower than competitors' plate profile.

## Bibliography

- Askew MJ, Mow VC, Wirth CR, Cambell CJ. Analysis of the intraosseus stress field due to compression plating. *J Biomech.* 1975; 8:203-12.
- Baggott DG, Goodship AE, Lanyon LE. A quantitative assessment of compression plate fixation in vivo: An experimental study using the sheep radius. *J Biomech.* 1981;14 (10):701-711.
- Beaupre GS, Carter DR, Dueland RT, Caler WE, Spengler DM. A biomechanical assessment of plate fixation, with insufficient bony support. *J Orthop Res.* 1988; 6:721-729.
- Carter DR, Vasu R, Spengler DM, Dueland RT. Stress fields in the unplated and plated canine femur calculated from in vivo strain measurements. *J Biomech.* 1981;14:63-70.
- Cordey J, Florin P, Klaue K, Perren SM. Compression achieved with the dynamic compression plate: effects of the inclined sloping cylinder and inclination of the screw. In: Uthoff HK, ed, *Current Concepts of Internal Fixation of Fractures.* Berlin: Springer-Verlag; 1980:192-200.
- Hayes WC, Perren SM. Flexural rigidity of compression plate fixation of fractures. *Proceedings of the 2nd Nordic Meeting on Medical and Biological Engineering.* 1971;242.
- Johansson JE, McBroom R, Barrington TW, Hunter GA. Fracture of the ipsilateral femur in patients with total hip replacement. *J Bone Joint Surg.* 1981;63-A (9):1435-1441.
- Klaue K, Perren SM, Kowalski M. Internal fixation with a self-compressing plate and screw: improvements of the plate hole and screw design, I. Mechanical investigation *J Orthop Trauma.* 1991;5:280-288.
- Lanyon LE, Hampson WGJ, Goodship AE, Shah JS. Bone deformation recorded from in vivo strain gages attached to the human tibial shaft. *Acta Orthop Scand.* 1975; 46:256-268.
- Nunamaker DM, Perren SM. A radiological histological analysis of fracture healing using prebending of compression plate. *Clin Orthop.* 1979;138:167-174.
- Serocki JH, Chandler RW, Coor LD. Treatment of fractures about hip prostheses with compression plating. *J of Arthroplasty.* 1992; (2):129-135.
- Zenni EJ, Pomeroy DL, Caudle, RJ. Ogden plate and other fixations for fractures complicating femoral endoprotheses. *Clin Orthop.* 1988 231:83-90.
- ISO 9585, Implants Surgery: Determination of bending strength and stiffness of bone plate.
- ASTM F 382-95, Standard Test Method for: Static bending properties of metallic bone plates.

