

## Kinectiv Technology Neck-Heads and Neck-Stems

Indicates "Combination is approved by Zimmer, Inc."

### Stainless Steel Femoral Heads, 12/14 Taper<sup>1</sup>

PROTASUL-S30 Femoral Head, +0 Length

Straight Necks 00-7848-002/044-00	Anteverted Necks 00-7848-002/044-01	Retroverted Necks 00-7848-002/044-01
✓	✓	✓

### CoCr Femoral Heads, 12/14 Taper<sup>1</sup>

DURASUL CoCr Femoral Head, +0 Length

METASUL Femoral Head +0 Length

METASUL LDH (Large Diameter Head)

METASUL LDH Head Adapter, +0 Length (12/14 taper)

VerSys CoCr Femoral Head, +0 Length

✓	✓	✓
✓	✓	✓
✓	✓	✓
✓	✓	✓
✓	✓	✓

### Ceramic Femoral Heads, 12/14 Taper<sup>1</sup>

BIOLOX-delta Femoral Head, +0 Length

BIOLOX-forte Femoral Head, +0 Length

BIOLOX-option Femoral Head, +0 Length

CERASUL Femoral Head System, +0 Length

SULOX Femoral Head, +0 Length<sup>2</sup>

Alumina 12/14 Ceramic Heads, +0 Length

✓	✓	✓
✓	✓	✓
✓	✓	✓
✓	✓	✓
✓	✓ <sup>2</sup>	✓ <sup>2</sup>
✓	✓	✓

### Femoral Stems<sup>1,2</sup>

M/L Taper Femoral Stem with Kinectiv Technology

CLS Brevius Stem with Kinectiv Technology

✓	✓	✓
✓	✓	✓

**1. Any Zimmer hip prosthesis with Kinectiv modular neck technology is designed to be used with a +0 head only. Use of longer heads could lead to fatigue failure of the stem or neck and use of shorter heads could lead to reduced range of motion and component impingement.**

**2. Due to component impingement, do not use SULOX ceramic heads with the following anteverted or retroverted modular necks:**

**00-7848-002-01 Modular Neck K1**

**00-7848-003-01 Modular Neck S1**

**00-7848-020-01 Modular Neck K2**

**00-7848-030-01 Modular Neck S2**